

Petri Dish

Hycan
Plastics



- Sterile Petri Dishes are produced from FDA certified virgin polystyrene and by ISO 13485 certified manufacturer
- CE marked in accordance with European Directive 98/79/EC
- Have three venting ribs into the underside of the lid to prevent condensation
- Available in one and two compartments
- Sterile by Ethylene Oxide (ETO)

Product code	Item Description	Vents	Base Dia.mm	Overall H.mm	Packaging
HCK-1004-PET	90x15 mm	Triple	89	16.2	500 pairs
HCK-1003-PET	90x15 mm, 2 compartments	Triple	89	16.2	500 pairs